

COVID-19 TROLLEY SET-UP

AIRWAY

- Emergency intubation checklist
- Team brief for intubating COVID-19 patient
- COVID-19 ventilator set-up guide
- Alternative laryngoscope
- ETT size 7 + 8 (uncut in case of proning)
- 10ml syringe
- Bougie
- Capnography
- HME filter (Covidien) DAR large
- Water's circuit - set up as per diagram
- Facemask
- Aquagel
- Chest drain clamp
- Tube tie
- 2 x Sputum traps

AIRWAY EQUIPMENT IN BOX

(Box can be cleaned and reused)

- iGel size 4 + 5
- Guedel airways
- Size 6 ETT
- Scalpel - No. 10 blade
- Magill forceps



LINES

- 2 x floswitch A-lines
- 2 x Chloroprep
- Gauze
- 3 x venflon dressings
- 2 ml syringe
- Orange needle
- 4 x 3-way taps
- 18G cannula
- Central line insertion pack (ensure 20cm quin-lumen)
- C-arm cover (for USS machine)
- Aquagel
- A-line transducer
- 500ml N. Saline for flush
- Giving set
- Pressure bag
- 3 x extension lines

OTHER:

- NG feeding tube + sticker
- Catheter
- Timer/stopwatch
- Green scissors
- 4 x viral swabs
- Blood culture pack
- Baby monitor
- Micropore tape

FURTHER ITEMS (RUNNER TO COLLECT)

- C-MAC + difficult airway trolley
- 2 x 10ml 1% lidocaine
- 2 x Sterile gloves (check size)
- USS machine
- Consider vascath
- Ambubag outside room

Please complete this list for every trolley made and leave on the trolley.